Gebruikershandleiding
Manuel d'utilisation

Smart Chair

! HORIZON

Geachte reiziger,

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LITHIUM-ION BATTERY

Unit Type: HQ24-12
 Nominal Voltage: 24v
 Rated Capacity: 288wh/12Ah
 Rated Charge Current: 2A

Important Note

Read the following usage and handling instructions for the lithium-ion batteries contained in the wheelchair. Our company will not be responsible for any injury or consequences resulting from failure to follow these instructions.

***DO NOT** use unauthorized chargers, as they may result fires or other potential risks.

***DO NOT** place battery near fire, heat, water or heavy moisture.

***DO NOT** disassemble or apply heavy force to the batteries.

***AVOID** short circuited the batteries by crossing the positive and negative terminals.

***DO NOT** mix different battery models or use batteries from other manufacturers.

***DO NOT** charge batteries for more than 8 hours or in high temperatures or severe cold

1

Important Note

environments.

*Keep batteries out the reach of children.

***STOP USING** the wheelchair and/or charging the batteries if the battery gets hot, emits an unusual odor or the casing gets compromised.

*If battery leaks due to mishandling, place in a safe area away from heat to prevent the batteries from catching fire.

*If the battery leaks fluid and makes contact with your eyes, wash immediately with water and avoid rubbing it into your eyes. Seek immediate medical attention to avoid serious injury.

*Use and storage - If discharging, charging and storing conditions are not followed, the battery performance and battery lifespan will be reduced and the battery may overheat leading to fire or explosion.

*Only use the battery charger supplied with the wheelchair to avoid potential risks.

2

Read the following usage and handling instructions for the lithium-ion batteries contained in the wheelchair. Our company will not be responsible for any injury or consequences resulting from failure to follow these instructions.

***DO NOT** pierce battery with a sharp object such as needles.

***DO NOT** hit, step on, throw, or drop the battery, as this may cause strong shocks.

***DO NOT** use a battery with serious scars or deformation.

*While the battery pack is being charged, used, or stored, keep it away from objects or materials that may have static electric charges.

*When not using the battery for an extended period, remove it from the equipment and store it in a place with low humidity and temperature.

*The battery can be used within the following temperature ranges.
DO NOT exceed these ranges:
 - Charge temperature range: **0°C to 45°C**
 - Discharge temperature range: **-20°C to 60°C**.

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LITHIUM-ION BATTERY

Unit Type: HQ24-12
Nominal Voltage: 24v
Rated Capacity: 288wh/12Ah
Rated Charge Current: 2A

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- ***DO NOT** use unauthorized chargers, as they may result fires or other potential risks.
- ***DO NOT** place battery near fire, heat, water or heavy moisture.
- ***DO NOT** disassemble or apply heavy force to the batteries.
- ***AVOID** short circuited the batteries by crossing the positive and negative terminals.
- ***DO NOT** mix different battery models or use batteries from other manufacturers.
- ***DO NOT** charge batteries for more than 8 hours or in high temperatures or severe cold

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*The battery can be used within the following temperature ranges.
DO NOT exceed these ranges:
- Charge temperature range: **0°C to 45°C**
- Discharge temperature range: **-20°C to 60°C**.

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

Joystick Malfunction Codes

Code	Information
E1	Upper and lower control communication failure
E2	Push handle to power on
E3	Overvoltage
E4	Upper control joystick abnormality
L1	Left motor electromagnet abnormality
L2	Left motor phase line abnormality
L3	Left motor Hall sensor abnormality
L4	Left motor phase line short circuit
J1	Right motor electromagnet abnormality
J2	Right motor phase line abnormality
J3	Right motor Hall sensor abnormality
J4	Right motor phase line short circuit

The image shows a joystick controller with a digital display. The display shows '00'. Below the main image are three inset images showing different display codes: 'L2', 'J2', and 'J4'. Each inset image also shows the 'SOFT', 'NORMAL', and 'SPORT' mode indicators.

IP24	waterproof and dustproof resist
	Body Floating Equipment
	Limit of Humidity
	Waste Electrical and Electronic Equipment
	Medical Device
	Consult Instructions for Use
	Do not use the device after it is damaged
	General Warning Sign
	Manufacturer
	Unique Device Information

manufacturer

SUZHOU KD INTELLIGENT DEVICE CO., LTD.

NO.36, GuGang Road, ChengXiang Town TaiCang City Jiangsu, CN
215400

EMAIL: Yiwenz@yahoo.com

Material Safety Data Sheet (MSDS) for passengerships and airlines.

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AB-AB006-a

任何未经授权更改、伪造或篡改本文件均属违法行为，将会受法律制裁。如有疑问，可扫二维码下载报告验证或者电话咨询。
Any unauthorized modification, forgery, or falsification of this document constitutes a violation of law and is subject to legal penalties. If you have any questions, you may scan the QR code to download the report for verification or contact us via email.

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3. Hazards Summarizing 危险概述**Danger sort 危险类别:** N/A**Routes of entry 进入途径:**

1. Eyes and Skin – When leaking, the electrolyte solution contained in the battery irritates to ocular tissues and the skin.
眼睛和皮肤 – 当电池泄漏时, 电池内部的电解液会刺激眼膜和皮肤, 甚至有疼痛感。
2. Inhalation – Respiratory (and eye) irritation may occur if fumes are released due heat or an abundance of leaking batteries.
吸入 – 电池大量泄漏产生热量导致冒烟, 吸入会刺激呼吸系统。
3. Ingestion – The ingestion of the battery can be harmful. Content of open battery can cause serious chemical burns of mouth, esophagus and gastrointestinal tract.
吞食 – 吞食电池对身体有很大伤害。电池里含的物质会引起嘴、食道和胃肠道化学灼伤。

Health harm 健康损害:

Exposure to leaking electrolyte from ruptured or leaking battery can cause 电池破裂导致电解液外漏会导致以下伤害:

1. Inhalation – Burns and irritation of the respiratory system, coughing, wheezing, and shortness of breath.
吸入 – 灼伤或刺激呼吸系统, 可能会产生咳嗽、喘息和呼吸浅短等现象。
2. Eyes – Redness, tearing, burns. The electrolyte is corrosive to all ocular

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AB-AB-006-a

任何未经授权更改、伪造或篡改本文件均属违法行为, 将会受法律制裁。若有疑问, 可扫二维码下载报告验证或者邮件咨询。
Any unauthorized modification, forgery, or falsification of this document constitutes a violation of law and is subject to legal penalties. If you have any questions, you may scan the QR code to download the report for verification or contact us via email.

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3. **Skin**—The electrolyte is corrosive and causes skin irritation and burns.
皮肤—电解液有腐蚀性, 会刺激皮肤甚至灼伤皮肤。4. **Ingestion**—The electrolyte solution causes tissue damage to throat and gastrointestinal track.

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4. First Aid Measures 急救措施

Skin contact: Not anticipated. If the battery is leaking and the contained material contacts the skin, flush with copious amounts of clear water for at least 15 minutes.**Eye contact:** Not anticipated. If the battery is leaking and the contained material contacts eyes, flush with copious amounts of clear water for at least 15 minutes. Get medical attention at once.**Inhalation:** Not anticipated. If the battery is leaking, remove to fresh air. If irritation persists, consult a physician.**Ingestion:** Not anticipated. If the battery is leaking and the contained material is ingested, rinse mouth and surrounding area with clear water at once. Consult a physician immediately for treatment.**皮肤接触:** 没有事先预料的, 如果电池漏液接触到皮肤上, 立即用大量的清水冲洗至少 15 分钟。**眼睛接触:** 没有事先预料的, 如果电池漏液接触到眼睛上, 立即用大量的清水冲洗至少 15 分钟, 并立即就医。**吸入:** 没有事先预料的电池泄漏, 转移到空气新鲜的地方, 如果刺激性还存在, 请咨询医生。**吞食:** 没有事先预料的, 如果电池漏液并且吞食了电池原料, 立即用清水冲洗嘴部及周围部位, 并就医治疗。

5. Fire Fighting Measures 消防措施

Unusual Fire and Explosion Hazards: Battery may explode or leak potentially hazardous vapors subject to: exposed to excessive heat (above the maximum rated temperature as specified by the manufacturer) or fire, over-charged, short circuit, punctured and crushed.**Hazardous Combustion Products:** Fire, excessive heat, or over voltage conditions may produce hazardous decomposition products. Damaged batteries can result in rapid heating and the release of flammable vapors.

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AB-AB-006-a

任何未经授权更改、伪造或篡改本文件均属违法行为, 将会受法律制裁。如有疑问, 可扫二维码下载报告验证或者邮件咨询。
Any unauthorized modification, forgery, or falsification of this document constitutes a violation of law and is subject to legal penalties. If you have any questions, you may scan the QR code to download the report for verification or contact us via email.

Hotline: 400-003-0500 web: www.anbotek.com.cn E-mail: service@anbotek.com

Anbotek

Product Safety

Report No.: 1811C50257312704
报告编号: 1811C50257312704Page 5 of 11
第 5 页 共 11 页

Extinguishing Media: Dry chemical type extinguishers are the most effective means to extinguish a battery fire. A CO₂ extinguisher will also work effectively.

Fire Fighting Procedures: Use a positive pressure self-contained breathing apparatus if batteries are involved in a fire. Full protective clothing is necessary. During water application, caution is advised as burning pieces of flammable particles may be ejected from the fire.

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6. Accidental Release Measures 意外泄漏措施

The material contained within the battery would only be released under abusive conditions. In the event of battery rupture and leakage, collect all the released materials that are not hot or burning in an appropriate waste disposal container while wearing proper protective clothing and ventilate the area. Placed in approved container and disposed according to the local regulations.

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7. Handling and Storage 操作和贮存

Handling 操作:

1. Batteries are designed to be recharged. However, improperly charging a battery may cause the battery to flame. When charging the battery, use dedicated chargers and follow the specified conditions.

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2. Never disassemble or modify a battery. 不拆解电池。
3. Do not immerse, throw, and wet a battery in water. 不浸没、投掷和用水弄湿电池。
4. Should a battery unintentionally be crushed, thus releasing its contents, rubber gloves must be used to handle all battery components. Avoid the inhalation of any vapors that may be emitted. 如果电池被无意挤压而导致内部物质释放, 必须戴上橡胶手套处理所有的电池成分, 避免吸入释放的任何气体。
5. Short circuit causes heating. In addition, short circuit reduces the life of the battery and can lead to ignition of surrounding materials. Physical contact with short-circuited battery can cause skin burn. 短路会引起电池过热。此外, 短路会使电池寿命大大减少, 甚至会导致周围材料着火。身体接触短路的电池会导致皮肤灼伤。

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AB-AB-006-a

任何未经授权更改、伪造或篡改本文件均属违法行为, 将会受法律制裁。如有疑问, 可扫二维码下载报告验证或者邮件咨询
Any unauthorized modification, forgery, or falsification of this document constitutes a violation of law and is subject to legal penalties. If you have any questions, you may scan the QR code to download the report for verification or contact us via email

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- Avoid reversing the battery polarity, which can cause the battery to be damaged or flame. 避免颠倒电池极性, 可能会引起电池损坏或者燃烧。
- In the event of skin or eye exposure to the electrolyte, refer to Section 4, First Aid Measures. 如果皮肤或者眼睛接触到电解液, 参考第四项并立即采取急救措施。

Storage 贮存:

- Batteries should be separated from other materials and stored in a noncombustible, well ventilated, sprinkler-protected structure with sufficient clearance between walls and battery stacks. Do not place batteries near heating equipment, nor expose to direct sunlight for long periods. 电池应该和其他材料分开并且贮存于通风且不易燃烧的地方。自动灭火装置应与墙和电池组保持足够的间隙。不要把电池靠近加热装置, 或者直接长时间的暴露于阳光直射的区域。
- Do not store batteries above 35°C or below -20°C. Store batteries in a cool (about 20±5°C) in a long time, dry and ventilated area that is subject to little temperature change. Elevated temperatures can result in reduced battery cycle life. Battery exposure to temperatures in excess of 60°C will result in the battery venting flammable liquid and gases. 不要在35°C以上和-20°C以下的环境贮存电池。电池应该贮存在干燥的、通风良好的阴凉区域(大约20±5°C)。升高温度会导致电池循环寿命减少。电池暴露于60°C以上的温度可能会导致电池泄漏可燃性液体和气体。
- Keep batteries in original package until use and do not jumble them. 保持电池最原始的包装直到使用时, 不要把电池弄混乱。

8. Exposure Controls/Personal Protection 暴露控制/自我防护

Engineering Controls: Keep away from heat and open flame.

Ventilation: Not necessary under conditions of normal use. In case of abuse, use adequate mechanical ventilation (local exhaust) for the battery that vent gas or fumes.

Respiratory Protection: Not necessary under conditions of normal use. If battery is burning, leave the area immediately. During fire fighting fireman should use self-contained breathing, full-face respiratory equipment. Fires may be fought but only from safe fire fighting distance, evacuate all persons from the area of fire immediately.

Eye Protection: Not necessary under conditions of normal use. Use safety glasses with side shields if handling a leaking or ruptured battery.

Body Protection: Not necessary under conditions of normal use. Use rubber apron and protective working in case of handling a leaking or ruptured battery.

Protective Gloves: Not necessary under conditions of normal use. Use chemical resistant rubber gloves if handling a leaking or ruptured battery.

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AB-AB-006-a

任何未经授权更改、伪造或篡改本文件均属违法行为, 将会受法律制裁。若有疑问, 可扫二维码下载报告验证或者邮件咨询。Any unauthorized modification, forgery, or falsification of this document constitutes a violation of law and is subject to legal penalties. If you have any questions, you may scan the QR code to download the report for verification or contact us via email.

Hotline: 400-003-0590 web: www.anbotek.com.cn E-mail: service@anbotek.com

Anbotek

Product Safety

Report No.: 1811C50257312704
报告编号: 1811C50257312704Page 7 of 11
第 7 页 共 11 页

Others: Use good chemical hygiene practice. Wash hands thoroughly after cleaning-up a battery spill caused by leaking battery. No eating, drinking, or smoking in battery storage area.

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10. Stability and Reactivity 稳定性和反应活性

Stability: Stable

Conditions to Avoid: Do not heat, throw into fire, disassemble, short circuit, immerse in water or overcharge, etc.

Incompatibility: None during normal operation. Avoid exposure heat, open flame and corrosives.

Hazardous Polymerization: Will not occur.

Hazardous Decomposition Products: The battery may release irritative gas once the electrolyte leakage.

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The battery does not elicit toxicological properties during routine handling and use. If the battery is opened through misuse or damage, discard immediately. Internal components of cell are irritant and sensitization.

Irritancy: The electrolytes contained in this battery can irritate eyes with any contact. Prolonged contact with the skin or mucous membranes may cause irritation.

Sensitization: No information is available.

Teratogenicity: No information is available.

Carcinogenicity: No information is available.

Mutagenicity: No information is available.

Reproductive toxicity: No information is available.

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12. Ecological Information 生态信息

1. When properly used and disposed, the battery does not present environmental hazard.

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2. The battery does not contain mercury, cadmium, or lead.

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3. Do not let internal components enter marine environment. Avoid releasing to water ways, wastewater or ground water.

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13. Disposal Considerations 废弃处理

1. Disposal of the battery should be performed by permitted, professional disposal firms knowledgeable in Federal, State or Local requirements of hazardous waste treatment and hazardous waste transportation. 处理电池要有许可, 在联邦、国家或者当地危物质处理部门和危物质运输部门要求的专业处理知识。
2. The battery should be completely discharged prior to disposal and/or the

Anbotek

Product Safety

Report No.: 1811C50257312704
报告编号: 1811C50257312704Page 9 of 11
第 9 页 共 11 页

terminals taped or capped to prevent short circuit. When completely discharged it is not considered hazardous. 处理电池之前要完全放电或者把电池末端用胶带粘上防止短路。完全放电的电池被认为是没有危害的。

3. The battery contains recyclable materials. Recycling options available in your local area should be considered when disposing of this product, through licensed waste Carrier.
电池包含可循环利用的材料。在当地回收利用这些处理掉的产品时,要取得废弃物处理的授权。

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AB-AB-006-a

任何未经授权更改、伪造或篡改本文件均属违法行为,将会受法律制裁。如有疑问,可扫二维码下载报告验证或者邮件咨询。
Any unauthorized modification, forgery, or falsification of this document constitutes a violation of law and is subject to legal penalties. If you have any questions, you may scan the QR code to download the report for verification or contact us via email
Hotline: 400-002-0500 web: www.anbotek.com.cn E-mail: service@anbotek.com

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ADR	European Agreement concerning the International Carriage of Dangerous Goods by Road 《国际危险货物道路运输欧洲协定》(ADR)相关规定
RID	Regulations concerning the International Carriage of Dangerous Goods by Rail 《国际危险货物铁路运输欧洲协定》(RID)相关规定

The dangerous goods regulations require that each battery design be subject to tests contained in Section 38.3 of the UN Manual of Tests and Criteria prior to being offered for transport.
危险品法规规定, 运输前, 每一个电池设计须通过联合国试验和标准手册38.3节所载的测试。

15. Regulatory Information 监管信息

- 《Dangerous Goods Regulations》
《危险物品规则》
- 《Recommendation on the Transport of Dangerous Goods Model Regulations》
《对危险货物运输的有关规定的建议》
- 《International Maritime Dangerous Goods》
《国际海运危险货物规则》
- 《Technical Instructions for the Safe Transport of Dangerous Goods》
《危险品安全运输技术指令》
- 《Classification and code of dangerous Goods》
《危险货物分类和品名编号》
- 《Consumer Product Safety Act》(CPSA)
《消费产品安全法》
- 《Federal Environmental Pollution Control Act》(FEPCA)
《联邦环境污染控制法》
- 《Resource Conservation and Recovery Act》(RCRA)
《资源保护及恢复法案》
- 《European Agreement concerning the International Carriage of Dangerous》
《国际危险货物道路运输欧洲协定》
- 《Regulations concerning the International Carriage of Dangerous》
《国际危险货物铁路运输欧洲协定》

In according with all Federal, State and local laws.
根据所有联邦、州和地方法律。

16. Other Information 其他信息

The information above is believed to be accurate and represents the best information currently available to us. However, this document makes no warranty of merchantability or any other warranty, express or implied, with respect to such information, and we assume no liability resulting from its use. Users should make their own investigations to determine the suitability of the information for their particular purposes. Although reasonable precautions have been taken in the

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preparation of the data contained herein, it is offered solely for your information, consideration and investigation. This material safety data sheet provides guidelines for the safe handling and use of this product; it does not and cannot advise on all possible situations, therefore, your specific use of this product should be evaluated to determine if additional precautions are required.

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* The MSDS is prepared based on the information provided by client. The contents and formats of this MSDS are revised as per client's request.

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AB-AB-006-a

任何未經授權更改、偽造或篡改本文件均屬違法行為, 將會受法律制裁。若有疑問, 可扫二维码下載報告驗證或郵件諮詢。
Any unauthorized modification, forgery, or falsification of this document constitutes a violation of law and is subject to legal penalties. If you have any questions, you may scan the QR code to download the report for verification or contact us via email.

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Shenzhen Anbotek Compliance Laboratory Limited
AB-AB-004-a